 ITSETESTI CRP-TASON MÄÄRITTÄMISEEN SEMIKVANTITATIIVISESTI  SJÄLVTEST FÖR ATT MÄTA SEMI-KVANTITATIVT CRP-HALTEN
 SELF-TEST FOR CRP CONCENTRATION DETERMINATION (SEMI-QUANTITATIVE)
 СДЕЛАЙТЕ ПОЛУКОЛИЧЕСТВЕННЫЙ АНАЛИЗ СРБ В КРОВИ

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INSTRUCTIONS OF USE

Thank you for choosing the Medigoo.com health quick test. This CRP (C-Reactive Protein) Rapid Test is used for the detection of CRP in blood to diagnose bacterial infectious disease and inflammatory disorders.

CRP is a non-specific marker which is mainly produced by the liver. In healthy patients, CRP concentration is lower than 8 mg/L while the concentration level can go higher than 100 mg/L in case of severe infection or during inflammatory process. Intermediate levels, within 8 and 100 mg/L, are concomitant with more or less mildly viral or bacterial infections that can be easily overcome by appropriate treatment ordered by your doctor. With respect to this test and CRP concentration, different lines will appear to the test cassette allowing the semi-quantitative interpretation of the test result.

PRODUCT CONTENTS

The box contains all necessary items to perform a test:

- 1 sealed aluminium pouch containing
- 1 test device, 1 plastic pipette, 1 desiccant bag
- 2 sterile lancet for blood sampling
- 1 processing bottle containing 3 ml of diluent
- 1 instruction leaflet

Open the protective aluminium pouch only when you are ready to perform the test. The desiccant bag should not be used.

Not provided necessary materials: clock for timing use.

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PRECAUTIONS

1. This test is exclusively intended for in vitro diagnostic only. External use only. DO NOT SWALLOW.
2. Carefully read the instructions before performing the test. The test is only interpretable if the instructions are carefully followed. Follow strictly the indicated blood and diluent quantities, as well as reading times.
3. Store between +4° C and +30° C. Do not freeze.
4. Do not use after the expiry date printed on the label and on the protective pouch or if this pouch is damaged.
5. Do not reuse test.
6. Keep out of reach of children.
7. After use, all the components can be discarded with other household garbage.

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QUESTIONS AND ANSWERS

How does the CRP-test work?

If present in the sample, CRP will interact with the anti CRP antibodies as well as with the CRP antigen placed in the test cassette producing different lines allowing a semi-quantitative measurement of the CRP concentration in the blood sample.

When can this test be performed?

CRP-test should be performed in the case of acute infection symptoms such as feverishness, fever, headaches or weakness. The test can be performed at any time of the day.

Can the result be incorrect?

The results are accurate as far as the instructions are carefully followed. Nevertheless, the result can be incorrect if CRP-test gets wet before test performing or if the quantity of diluted blood dispensed in the sample well is not sufficient. The plastic pipette provided in the box allows making sure the collected blood volume is correct and is correctly diluted into the diluent dropper vial.

How to interpret the test if the colour and the intensity of the lines are different?

The colour intensity of the lines is important to interpret the test result. So please refer to the pictures and instructions in paragraph "Result interpretation".

What is the line that appears under the mark C (Control) for?

When this line appears, it only means that the test was performing well.

If I read the result after 10 minutes, will the result be reliable?

No. The result should be read within 10 minutes. The result is reliable up to 10 minutes but should be better read at 5 minutes after adding the sample.

What do I have to do if the result is positive?

The result is positive when the CRP is higher than 8 mg/L, please refer to the result interpretation and inform your doctor as soon as possible as an antibiotic treatment might be urgently needed.

What do I have to do if the result is negative?

The result is negative when the CRP is lower than 8 mg/L, it means that there is neither bacterial nor viral infection on the way. However if the symptoms persist, it is recommended to consult your doctor.

What is the accuracy of CRP-test ?

The CRP-test is accurate and is used for more than 10 years by professionals in the field (hospitals, laboratories...).

Evaluation reports show that an overall agreement of at least 98% with reference method. Although this test is reliable, false positive or false negative results can be obtained.

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🇬🇧 Wash your hands thoroughly. Use soap and warm water. Dry your hand with a clean towel

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🇬🇧 Open the foil pouch (from the notch) and only get out the device and the pipette. Dispose of the desiccant pouch.

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🇬🇧 Unscrew the processing bottle containing the diluent

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🇬🇧 Prepare the lancet. Twist off and discard the lancet's grey cap without touching the trigger button.

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🇬🇧 Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.

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🇬🇧 Strongly press the extremity of the lancet towards previously massaged ring or middle finger. Puncture the skin by pressing the trigger button.

The tip retracts automatically and safely after use.

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🇬🇧 Gently rub the hand from the wrist and palm down to the finger to form a rounded drop of blood over the puncture site.

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🇬🇧 Without pressing the pipette bulb, put it in contact with the drop. The blood migrates into the pipette through the capillarity to the line indicated on the pipette. You may massage again your finger to obtain more blood if the line is not reached. As far as possible avoid air bubbles. If the capillary effect is not working, you may suck in the blood into the pipette by pressing the pipette bulb.

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🇬🇧 Place the pipette into the processing bottle liquid and eject all the blood. Avoid bubbles. Close the cap. Turn the bottle upside down several times. Mix well, but no shaking.

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🇬🇧 Without opening the sample processing bottle, break the tip of the cap.

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🇬🇧 Drop 4 drops of processed specimen on the cassette. Read the result at 5 minutes. Do not interpret after 10 minutes.

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Only 2 lines appear (no line at the mark T). This result means that there is neither infection nor inflammation.

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Three colored lines appear in the window but the intensity of the line T mark is weaker than the intensity of the two other lines. This result means that an infection could be on the way. You should consult your doctor.

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Three colored lines appear in the window and the intensity of the three lines is more or less equal. This result means that a viral or bacterial infection is active or on its way. You must consult your doctor.

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Three colored lines appear in the window but the color intensity of the lines under T and C marks are stronger than the color intensity of the middle line. In case of very high CRP concentration, the middle line may even disappear. This result means a severe bacterial infection is on the way. You must consult your doctor urgently !!!

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If no line appears at the mark C, it is impossible to interpret the test. It must be considered as non-valid, even if a colored line is visible under T mark and/or the middle line appears. It is recommended to repeat the test with a new CRP Test.

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©Medigoo 2016-2018
Owen Mumford Ltd