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INSTRUCTIONS FOR USE

Thank you for choosing Medigoo.com personal health test kit. This PSA Prostate Specific Antigen Test is a rapid chromatographic immunoassay for the qualitative detection of Prostate Specific Antigen in blood.

In man, the prostate antigen (PSA) is secreted by the prostate. This gland, situated in the abdomen under the bladder, surrounds the initial part of the urethra and plays an important role in the semen production. The prostate antigen level gives information about the prostate physiological condition. Thus, levels above the norm may show pathology of the prostate (benign hypertrophy, prostatitis, cancer, etc.). The PSA assay could be performed at least once a year in men over 50 years. In case of prostate problem, early screening allows considerably increasing recovery possibilities.

PRODUCT CONTENTS

The box contains all necessary items to perform a test:

- 1 Sealed aluminium pouch containing:
 - 1 test device, 1 plastic pipette, 1 desiccant bag
 - 1 diluent bottle containing 1 mL of diluent
 - 2 sterile lancets for blood sampling
 - 1 instruction leaflet.

Open the protective aluminium pouch only when you are ready to perform the test. The desiccant bag should not be used.

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PRECAUTIONS

1. This test is exclusively intended for in vitro diagnostic only. External use. DO NOT SWALLOW.
2. Carefully read the instructions before performing the test. The test is only interpretable if the instructions are carefully followed. Follow strictly the indicated blood and diluent quantities, as well as reading times.
3. Store between +4° C and +30° C. Do not freeze.
4. Do not use after the expiry date printed on the label and on the protective pouch or if this pouch is damaged.
5. Do not reuse the test.
6. Keep out of reach of children.
7. After use, all the components can be discarded along with the household garbage disposal.

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QUESTIONS AND ANSWERS

How does PSA test work?

Prostate specific antigen is a protein secreted in the semen by the prostate. A part of this protein is found in blood. The blood PSA level allows an estimation of the prostate physiological condition. A value higher than the norm indicates prostate hypertrophy. This hypertrophy may be benign, may be due to prostatitis or may be cancerous in case of adenocarcinoma.

Prostate pathology screening is for men over 50 years and is performed up to 75 years. The PSA test allows one to determine a high blood level, above 4 ng/mL, of prostate specific antigen likely to show benign or malignant prostate hypertrophy.

When should this test be used?

The PSA test can be used at any time of the day. However, in the following cases, it is recommended to wait the indicated time before using the test, to avoid false positive results:

Waiting time

- Cycling/Ergometry 24 hours
- Ejaculation 24 hours
- Prostate massage 2-3 days
- Transrectal ultrasounds 2-3 days
- Cystoscopy 1 week
- Prostate transurethral resection / biopsy 4-6 weeks

Can the results be incorrect?

The results are accurate as long as the instructions are carefully followed. Nevertheless, the result can be incorrect if the PSA test gets wet before performing the test or if the quality of the blood put in the sample well is not sufficient. The plastic pipette provided in the box is designated to ensure the collected blood volume is correct.

How to interpret the test if the colour and the intensity of the lines are different?

The colour and intensity of the lines have no importance for result interpretation. The lines should only be uniform and full.

What is the line that appears under the C (Control) mark for?

When this line appears, it only means that the test has worked correctly.

If I read the result after 15 minutes, will the result be reliable?

No. The test should be read at 10 minutes of adding the diluent. The results are reliable up to 15 minutes.

What do I have to do if the result is positive?

If the result is positive, it means that the PSA level present in your blood is higher than the norm (4 ng/mL) and that you should consult a doctor to show him the test results. Then, the doctor will decide what you have to do.

What do I have to do if the result is negative?

If the result is negative, it means that the PSA level in your blood is below 4 ng/mL. For patients aged between 50 and 75 years or having relatives suffering from prostate cancer, it is recommended to perform a test every year.

What is the accuracy of PSA test?

The PSA test is accurate and has been used for more than 10 years by professionals (hospitals, laboratories) in the field. Evaluation reports show an overall agreement of at least 88% with reference methods. Although this test is reliable, false positive or false negative results could be obtained.

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Thank you for choosing MediGoo.com !

1

Wash your hands thoroughly with soap and warm water. Dry your hands with a clean towel.

Wash your hands thoroughly. Use soap and warm water. Dry your hands with a clean towel.

Wash your hands thoroughly. Use soap and warm water. Dry your hands with a clean towel.

Wash your hands thoroughly. Use soap and warm water. Dry your hands with a clean towel.

2

Open the foil pouch (from the notch) and only use the test cassette and the pipette. Dispose of the desiccant pouch.

Open the foil pouch (from the notch) and only use the test cassette and the pipette. Dispose of the desiccant pouch.

Open the foil pouch (from the notch) and only use the test cassette and the pipette. Dispose of the desiccant pouch.

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3

Prepare the lancet. Twist off and discard the lancet's grey cap without touching the trigger button.

Prepare the lancet. Twist off and discard the lancet's grey cap without touching the trigger button.

Prepare the lancet. Twist off and discard the lancet's grey cap without touching the trigger button.

Prepare the lancet. Twist off and discard the lancet's grey cap without touching the trigger button.

4

Massage the hand without touching the future puncture site by rubbing down the hand towards the fingertip of the middle or ring finger. Choose the finger you use the least. You can clean the fingertip by a wipe damped in disinfectant or alcohol.

Massage the hand without touching the future puncture site by rubbing down the hand towards the fingertip of the middle or ring finger. Choose the finger you use the least. You can clean the fingertip by a wipe damped in disinfectant or alcohol.

Massage the hand without touching the future puncture site by rubbing down the hand towards the fingertip of the middle or ring finger. Choose the finger you use the least. You can clean the fingertip by a wipe damped in disinfectant or alcohol.

the fingertip of the middle or ring finger. Choose the finger you use the least. You can clean the fingertip by a wipe damped in disinfectant or alcohol.

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Strongly press the tip of the lancet towards previously chosen finger. Puncture the skin by pressing the trigger button. The tip retracts automatically and safely after use.

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Gently rub the hand from the wrist and palm down to the finger to form a rounded drop of blood over the puncture site.

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Without pressing the pipette bulb, put it in contact with the blood drop. The blood migrates into the pipette through the capillary effect to the line indicated on the pipette. You may massage again your finger to obtain more blood if the line is not reached. Take blood only upto the line on the pipette. Avoid air bubbles.

Without pressing the pipette bulb, put it in contact with the blood drop. The blood migrates into the pipette through the capillary effect to the line indicated on the pipette. You may massage again your finger to obtain more blood if the line is not reached. Take blood only upto the line on the pipette. Avoid air bubbles.

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Put all the collected blood from the pipette into the round sample well (s) of the cassette, by pressing on the pipette bulb.

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Wait 30-40 seconds for the blood to be totally absorbed into the sample well (s). Unscrew the blue cap of the plastic dropper bottle (while keeping the white cap tightly screwed in), turn the bottle upside down and add exactly 4 drops of diluent into the sample well (s) with an interval of 2-3 seconds between each drop.

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Read the results at 10 min. Do not interpret the results after 15 min.

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One colored line appears at the control region (C). This result means that the prostate specific antigen level is normal. A negative result cannot rule out the possibility of a prostate cancer. PSA level is below 4 ng/ml.

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Two distinct colored lines appear. One line should be in the control region (C) and another line should be in the test region (T). The intensity of the color in the test line region (T) will vary depending on the concentration of PSA present in the specimen. The result means that the prostate specific antigen level is higher than the norm and that you should consult a physician. PSA level is over 4 ng/ml.

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Test failed: One red line is shown in the Test area (T) and no line in the control area C. Alternatively, no lines are shown either at test area T or control area C. Repeat the test with a new testing device and with fresh blood.

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©Medigoo 2016-2018
FACTORY Owen Mumford Ltd