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INSTRUCTIONS FOR USE

Thank you for choosing Medigoo.com quick health test. This health test kit is intended to test the most common type of Anemia, iron deficiency Anemia. There are other Anemia types, which can not be detected by this test. Your ferritin level in the blood is tested. Ferritin is a protein that stores iron and iron deficiency can lead to anemia condition.

Anemia due to iron depletion is widely held in children and women of all ages but mainly in women who still have their period (at least 20% suffer iron deficiency). Main signs as paleness, feeling tired, headaches, faster heartbeat or shortness of breath during exercise appear gradually and could go unnoticed. Therefore it is important to determine if the available iron is sufficient for the body's needs.

Depletion in iron, which can happen during pregnancy, growth, in case of insufficient iron intake, inadequate absorption or blood loss (periods, abnormal bleedings, ulcers...etc) have tremendous effects on health.

The ferritin levels correspond to iron levels and low levels result to reduced production of hemoglobin. The red blood cells mainly consist of hemoglobin, which binds oxygen and is important in oxygen transport in the whole body. The lack of hemoglobin or red blood cells is anemia condition.

Caution: This test is not appropriated for patient suspected of or suffering from hemochromatosis.

PRODUCT CONTENTS

The box contains all necessary items to perform a test:

- 1 Sealed aluminium pouch containing:
 - 1 test device, 1 plastic pipette, 1 desiccant pouch
 - 1 diluent bottle containing 1 mL of diluent
 - 2 sterile lancet for blood sampling
 - 1 instruction leaflet.

Open the protective aluminium pouch only when you are ready to perform the test. The desiccant bag should not be used.

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QUESTIONS AND ANSWERS

How does Anemia test work?

Ferritin is a protein that binds and stores most of the iron in the body. Neither only a minor part of the ferritin is found in blood, it corresponds the amount of iron stored in the body. Positive result means that ferritin concentration in blood is lower than 20 ng/mL and a possible iron deficiency.

When can this test be performed?

Anemia test must be performed in case of symptoms like paleness, feeling tired, headaches, faster heartbeat or shortness of breath during exercise mainly when pregnant or in case of copious periods. The test can be used at any time of the day but must not be performed in case of diseases, acute inflammations or in case of spleen or liver injury. Positive results could be obtained even in case of no iron deficiency situation.

Can the result be incorrect?

The results are accurate as far as the instructions are carefully followed. Nevertheless, the result can be incorrect if the test gets wet before test performing or if the quantity of blood dispensed in the sample well is not sufficient. The plastic pipette provided in the box allows making sure the collected blood volume is correct.

How to interpret the test if the colour and the intensity of the lines are different?

The colour and intensity of the lines have no importance for result interpretation. The lines should only be homogeneous and clearly visible. The test should be considered as negative whatever the colour intensity of the test line.

What is the line that appears under the mark C (Control) for?

When this line appears, it only means that the test was performing well.

If I read the result after 15 minutes, will the result be reliable?

No. The result should be read at 10 minutes after adding the diluent. The result is reliable up to 15 minutes.

What do I have to do if the result is positive?

If the result is positive, it means that the ferritin level in blood is lower than the norm (20ng/mL) and that you should consult a doctor to show the test result. Then, the doctor will decide whether additional analysis should be performed.

What do I have to do if the result is negative?

If the result is negative, it means that the ferritin level is higher than 20 ng/mL and is within the norm. However, if the symptoms persist, it is recommended to consult a doctor.

What is the accuracy of Anemia test?

The test is accurate and has been used for more than 10 years by professionals in the field. Evaluation reports show an overall agreement of at least 98% with reference methods. Although this test is reliable, false positive or false negative results could be obtained.

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+ Wash your hands thoroughly. Use soap and warm water. Dry your hands with a clean towel.

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+ Open the foil pouch (from the notch) and only use the test cassette and the pipette. Dispose of the desiccant pouch.

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+ Prepare the lancet. Twist off and discard the lancet's grey cap without touching the trigger button.

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+ Massage the hand without touching the future puncture site by rubbing down the hand towards

the fingertip of the middle or ring finger. Choose the finger you use the least. You can clean the fingertip by a wipe damped in disinfectant or alcohol.

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+ Strongly press the tip of the lancet towards previously chosen finger. Puncture the skin by pressing the trigger button. The tip retracts automatically and safely after use.

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+ Gently rub the hand from the wrist and palm down to the finger to form a rounded drop of blood over the puncture site.

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+ Without pressing the pipette bulb, put it in contact with the blood drop. The blood migrates into the pipette through the capillary effect to the line indicated on the pipette. You may massage again your finger to obtain more blood if the line is not reached. Take blood only upto the line on the pipette. Avoid air bubbles.

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+ Put all the collected blood from the pipette into the round sample well (s) of the cassette, by pressing on the pipette bulb.

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+ Wait 30-40 seconds for the blood to be totally absorbed into the sample well (s). Unscrew the blue cap of the plastic dropper bottle (while keeping the white cap tightly screwed in), turn the bottle upside down and add exactly 4 drops of diluent into the sample well (s) with an interval of 2-3 seconds between each drop.

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+ Read the results at 10 min. Do not interpret the results after 15 min.

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➤ Two coloured lines appear in the window under the T (Test) and C (Control) marks. The intensity of the T line may be more transparent than the intensity of the C line. The result means that the ferritin concentration in blood is normal and that there is no potential iron deficiency.

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➤ Only one coloured line appears at the mark C (Control). This result means that the ferritin concentration in blood is too low. Iron reserves are insufficient. You should consult a doctor because it may be an iron deficiency anemia.

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➤ Test failed: One red line is shown in the Test area (T) and no line in the control area C. Alternatively, no lines shown either at test area T or control area C. Repeat the test with a new testing device and with fresh blood.

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©Medigoo 2016-2017
Owen Mumford Ltd